

# Total Intravenous Sets

Anaesthesia





## Why Choose Tivatek?

Absolute focus on creating a gold standard in Total Intravenous Anaesthesia, TIVA

Our motivation is to create an advanced range of TIVA products following guidelines from the Safe Anaesthesia Liaison Group (2009). We pride ourselves on exceptional manufacturing standards for a range of TIVA sets you can rely on. We aim to offer an all in one solution to make delivery of TIVA as safe as possible.



- Anti reflux valve on side arm for infusion of gravity fed IV solutions and anti siphon/anti reflux valves on drug lines (to prevent free-flow and back-flow) for the syringe pumps included.
- Microbore tubing used in all Tivatek
  TM kits ensures a low priming
  volume and reduces drug wastage.

- Anti-kink bridge as additional safety feature.
- High specification kink resistant tubing offers a reduced risk of occlusion.



### **Product Information: TT202E**

#### **TT202E**

### 2-way TIVA Infusion Set with high flow valve

Priming volume mL: 2.2 Carton size: 38 x 29 x 21cm Mandatory field Brand Name:

Units per box: 50 Carton quantity: 50 Europlaz - 2-way TIVA 2.0 m Infusion Set

Length: 2.0m Latex / DEHP?: No

Description:

2-way TIVA infusion set (Blue Caps)

Total length: 200cm

Lumen 1 - 'Side arm' with Non-Return valve for gravity fed infusion.

Lumen 2 - Drug Line with clamp and Pressure Activated Non-return / Anti-Syphon valve for Syringe-pump drug infusion.

Typical Application:

For single drug sedation or delivery of opioid analgesic agents.